

*Note: Modular TFF Units will vary in appearance

Performance Meets Flexibility: Modular TFF Benchtop System

Introducing Tezalon's groundbreaking Modular TFF Benchtop System. Designed to enhance your tangential flow filtration process, our system combines unparalleled flexibility with ease of use, empowering you to optimize your workflow.



Configurable, customizable, convenient

Are you focused on flexibility in TFF development and data generation? With Tezalon's Modular TFF, you can streamline your processes without getting bogged down in complicated setups. Why waste time learning intricate recipes when you can quickly analyze your data and pivot your TFF approach as needed?

Your System, Your Control, Your Choice

By partnering with best-in-class equipment from Levitronix, Quattroflow, Masterflex, Endress House, and Pendotech, we are able to connect the flexibility you require with the reliability of brands you trust.

Let's Talk Modular TFF!

Stages and scales, therapeutic modalities, evaluation of filter choices all require that you have an across-the-board non-dedicated system that permits you to configure what you need.

Let's be frank... with dozens of equipment vendors, endless and growing therapeutic modalities, and evolving process-critical parameters, don't be stuck with a proprietary control system, platform, sensor, or filter.

Invest in a system for development and pilot work that can be modified for future use on other therapeutic projects. Handle modality changes (Mab to CT to GT) and modern advances in process sciences with ease.

Choose a pilot scale vendor that understands configuration of dedicated skids for your mature process and helps you make the transition to production when you are ready.


Tablet controlled manual HMI, with the ability to click on valves and pumps to actuate process parameters. You control the process simply.

Experience unmatched control and versatility with Tezalon's Modular TFF Benchtop System.

Your TFF, your way—designed for seamless adaptability and optimal performance.

Unlock the full potential of your bioprocessing operations with Tezalon, where flexibility meets efficiency for unparalleled results.

Flexible components can be configured to any footprint and added/subtracted at will. Tube sets can be defined or users can create and modify to their needs. Valves, pumps, viscosity, concentration, UV, pressure, conductivity, flow, pH and PAT all possible.



For more information or product support:

+1.661.381.7950
info@tezalon.bio
www.tezalon.bio